


## Disposable 50-mm Sterile Petri Dishes, PS

### No need to presterilize

- Sterilized dishes are ready to use
- Use only what you need—come in convenient packs of 20
- Polystyrene dishes come with or without 47-mm absorbent cellulose pad that hold 1.8 to 2.0 mL of liquid media
- Ridges for easy stacking



Description	Dia x H (mm)	Qty/cs	Catalog number	Price/cs
Petri dish with pad	50 x 11	100	<a href="#">GY-14005-20</a>	
		500	<a href="#">GY-14005-28</a>	
Petri dish without pad	50 x 11	100	<a href="#">GY-14005-24</a>	
		500	<a href="#">GY-14005-26</a>	

[GY-14005-22](#) Presterilized absorbent pads. Box of 100 individual packages

## Petri Dishes, PFA

### Maintain purity of your culture

- Ideal for contamination-free biological applications
- Inert, nonstick, noncontaminating surface
- Lids not included
- Autoclavable



Capacity (mL)	Dia x H (mm)	Qty/pk	Catalog number	Price/pk
60	50 x 15	1	<a href="#">GY-06148-00</a>	
100	100 x 15	1	<a href="#">GY-06148-02</a>	

## Disposable Gridded Sterile Contact Plates, PS

### Identify and monitor cultures easily

- For determining surface microbial contamination
- Dishes have a raised agar bed for surface contact
- Stackable lids are flared for easy handling and automatic filling



Dia x H (mm)	Catalog number	Qty/cs	Price/cs
65 x 15	<a href="#">GY-14006-50</a>	500	

The grid consists of sixteen 10 x 10 mm squares. Disposable polystyrene (PS). Radiation sterilized.

## Petri Dishes, Glass

### Save money—autoclavable for repeated use

- Dishes withstand repeated wet or dry sterilization
- Easy sorting—tops marked in blue, bottoms marked in white



Dia x H (mm)	Catalog number	Qty/pk	Price/pk	Catalog number	Qty/cs	Price/cs
60 x 15	<a href="#">GY-34551-01</a>	12		<a href="#">GY-34551-21</a>	72	
100 x 10	<a href="#">GY-34551-02</a>	12		<a href="#">GY-34551-22</a>	72	
100 x 15	<a href="#">GY-34551-03</a>	12		<a href="#">GY-34551-23</a>	72	
100 x 20	<a href="#">GY-34551-04</a>	12		<a href="#">GY-34551-24</a>	72	
150 x 15	<a href="#">GY-34551-05</a>	12		<a href="#">GY-34551-25</a>	24	
150 x 20	<a href="#">GY-34551-06</a>	12		<a href="#">GY-34551-26</a>	24	

## Sterile Petri Dishes, PS

### Stacks won't tip over—convenient stacking ridges for better stability

- Use what you need—packaged in bags of 20 dishes
- Save time—already sterilized by gamma irradiation



Dia x H (mm)	Qty/pk	Catalog number	Price/pk
60 x 15	500	<a href="#">GY-06139-01</a>	
100 x 15	500	<a href="#">GY-06139-00</a>	
150 x 20	100	<a href="#">GY-06139-02</a>	

## Disposable Cell Spreaders

### Easy spreading with L-shaped design

- Won't cut into agar—upturned end to protect from puncture

Made of durable polystyrene, these sterile, disposable spreaders provide safe, convenient sample plating without flaming.



Description	Qty/pk	Qty/cs	Catalog number	Price/cs
Disposable cell spreaders	10	500	<a href="#">GY-14201-00</a>	
	1	500	<a href="#">GY-14201-10</a>	

**D**

**Dishes**

Dishes

**Crystallizing Dishes, Glass**

Reduced chipping—reinforced and fire-polished rims

- Withstand repeated sterilizations—wet or dry
- Autoclavable



Capacity (oz)	Approx dia x H (mm)	Catalog number	Qty/pk	Price/pk	Catalog number	Qty/pk	Price/cs
180 (6)	70 x 50	<a href="#">GY-34550-01</a>	6		<a href="#">GY-34550-21</a>	24	
180 (6)	80 x 40	<a href="#">GY-34550-02</a>	6		<a href="#">GY-34550-22</a>	24	
270 (9)	90 x 50	<a href="#">GY-34550-03</a>	6		<a href="#">GY-34550-23</a>	18	
325 (11)	100 x 50	<a href="#">GY-34550-04</a>	6		<a href="#">GY-34550-24</a>	18	
740 (25)	125 x 65	<a href="#">GY-34550-05</a>	4		<a href="#">GY-34550-25</a>	12	
1200 (30)	150 x 75	<a href="#">GY-34550-06</a>	4		<a href="#">GY-34550-26</a>	8	
1770 (60)	170 x 90	<a href="#">GY-34550-07</a>	2		<a href="#">GY-34550-27</a>	8	
2500 (85)	190 x 100	<a href="#">GY-34550-08</a>	2		<a href="#">GY-34550-28</a>	6	

**Standard-Form Evaporating Dishes, Porcelain**

Clean up easier with glazed interior

- Easy-pour spout
- Autoclavable



Capacity (mL)	Top outside dia (mm)	Height (mm)	Catalog number	Qty/pk	Price/pk	Catalog number	Qty/cs	Price/cs
35	60	24	<a href="#">GY-17812-00</a>	12		<a href="#">GY-17812-20</a>	48	
70	75	30	<a href="#">GY-17812-02</a>	12		<a href="#">GY-17812-22</a>	48	
80	80	30	<a href="#">GY-17812-04</a>	6		<a href="#">GY-17812-24</a>	24	
120	94	42	<a href="#">GY-17812-06</a>	6		<a href="#">GY-17812-26</a>	24	
150	100	42	<a href="#">GY-17812-07</a>	6		<a href="#">GY-17812-27</a>	18	
250	115	45	<a href="#">GY-17812-08</a>	6		<a href="#">GY-17812-28</a>	18	
385	145	48	<a href="#">GY-17812-09</a>	1		<a href="#">GY-17812-29</a>	12	
525	162	51	<a href="#">GY-17812-10</a>	1		<a href="#">GY-17812-30</a>	12	
765	185	54	<a href="#">GY-17812-12</a>	1		<a href="#">GY-17812-32</a>	8	
1285	215	75	<a href="#">GY-17812-14</a>	1		<a href="#">GY-17812-34</a>	6	
2100	265	80	<a href="#">GY-17812-16</a>	1		<a href="#">GY-17812-36</a>	4	
3250	305	95	<a href="#">GY-17812-18</a>	1		<a href="#">GY-17812-38</a>	3	

**Evaporating Dishes, Glass**

Reusable dishes save you money

- Flat-bottom evaporating dishes with a lip and spout and rugged beaded rims
- Autoclavable



Capacity (mL)	Top outside dia (mm)	Height (mm)	Catalog number	Price/pk	Catalog number	Qty/cs	Price/cs
140	80	45	—	—	<a href="#">GY-34586-20</a>	6	
170	90	50	—	—	<a href="#">GY-34586-22</a>	6	
290	105	55	<a href="#">GY-34586-04</a>		<a href="#">GY-34586-24</a>	6	

**Shallow-Form Evaporating Dishes, Porcelain**

Easy pouring spout for less mess

- Glazed interior and partial exterior
- Autoclavable



Capacity (mL)	Top outside dia (mm)	Height (mm)	Catalog number	Qty/pk	Price/pk	Catalog number	Qty/cs	Price/cs
30	70	16	<a href="#">GY-17814-00</a>	6		<a href="#">GY-17814-20</a>	24	
50	80	20	<a href="#">GY-17814-02</a>	6		<a href="#">GY-17814-22</a>	24	
75	90	22	<a href="#">GY-17814-04</a>	1		<a href="#">GY-17814-24</a>	18	
100	105	23	<a href="#">GY-17814-06</a>	1		<a href="#">GY-17814-26</a>	18	
195	120	34	<a href="#">GY-17814-08</a>	1		<a href="#">GY-17814-28</a>	12	
290	160	31	<a href="#">GY-17814-10</a>	1		<a href="#">GY-17814-30</a>	8	

**Petri Dish Racks**

Organize your dishes to save valuable bench space

- Easy viewing of incubating petri dishes
- Use with plastic dishes only

Capacity	Dimensions (L x W x H)	Cat. no.	Price/ea
56 dishes, 60 mm dia	10½" x 6½" x 5½" 267 x 165 x 140 mm	<a href="#">GY-06146-00</a>	
24 dishes, 100 mm dia	8½" x 4" x 8" 206 x 102 x 203 mm	<a href="#">GY-06146-10</a>	
60 dishes, 100 mm dia	13" x 7¾" x 7¾" 330 x 200 x 197 mm	<a href="#">GY-06146-30</a>	

