

Funnels

Powder Funnels, PP

Minimize powder bridging during transfer

- 60° wall incline
- Autoclavable



| Capacity mL (oz) | Top OD | Stem OD | Total height | Catalog number | Qty/pk | Price/pk |
|------------------|--------|---------|--------------|-----------------------------|--------|----------|
| 49 (1.7) | 65 mm | 15 mm | 71 mm | GY-06123-10 | 12 | |
| 138 (4.7) | 80 mm | 18 mm | 83 mm | GY-06123-20 | 12 | |
| 248 (8.4) | 100 mm | 21 mm | 97 mm | GY-06123-30 | 6 | |
| 748 (26.5) | 150 mm | 29 mm | 136 mm | GY-06123-40 | 4 | |

Vented Funnels, PP

No more air lock for easier funneling

- No messy cleanup—disposable for easier sampling
- Autoclavable



| Capacity mL (oz) | Top OD | Stem OD | Total height | Catalog number | Qty/pk | Price/pk |
|------------------|--------|---------|--------------|-----------------------------|--------|----------|
| 40 (1 1/3) | 58 mm | 6 mm | 70 mm | GY-06124-10 | 100 | |
| 50 (1 3/4) | 67 mm | 6 mm | 81 mm | GY-06124-20 | 100 | |

Large Industrial-Size Funnels, LDPE, PP, and HDPE

Quick, easy transfer of large quantities

- LDPE funnels have external ribbing for easy air release and a handle tab for easier lifting
- PP and HDPE funnels are translucent




| Capacity L (gal.) | Top OD | Stem OD | Total height | Catalog number | Qty/pk | Price/pk |
|---|---------|---------|--------------|-----------------------------|--------|----------|
| Low-density polyethylene (LDPE) | | | | | | |
| 1.5 (3/8) | 7 5/8" | 3/4" | 9 1/2" | GY-06121-10 | 1 | |
| 4 (1) | 11" | 7/8" | 11 1/2" | GY-06121-20 | | |
| Polypropylene (PP), autoclavable | | | | | | |
| 2 (1/2) | 8 3/16" | 1 1/16" | 7 7/8" | GY-06121-30 | 1 | |
| 4.5 (1 1/4) | 10 1/4" | 1 5/16" | 10 9/16" | GY-06121-32 | | |
| 10 (2 1/2) | 14" | 1 9/16" | 17 3/8" | GY-06121-34 | | |
| High-density polyethylene (HDPE) | | | | | | |
| 12 (3 1/4) | 15 3/8" | 1 7/8" | 14 3/8" | GY-06121-00 | 1 | |
| 18 (5) | 17 1/4" | 1 3/4" | 16 3/8" | GY-06121-36 | | |

Cole-Parmer® Utility Funnels, PP

Filter samples faster

- Eight internal ribs facilitate easier filtration
- Autoclavable



| Capacity mL (oz) | Top OD | Stem length | Qty/pk | Catalog number | Price/pk |
|-------------------|--------|-------------|--------|-----------------------------|----------|
| 5 (3/16) | 37 mm | 30 mm | 12 | GY-63100-26 | |
| 25 (3/4) | 52 mm | | | GY-63100-28 | |
| 40 (1 1/3) | 67 mm | | | GY-63100-30 | |
| One of each above | | | 3 | GY-63100-32 | |

Utility Funnels, PP and HDPE

Get easy, economical funneling

- Ideal for general use or powders
- PP versions are autoclavable



| Capacity mL (oz) | Top OD | Stem OD | Total height | Qty/pk | Catalog number | Price/pk |
|---|--------|---------|--------------|--------|-----------------------------|----------|
| Polypropylene (PP), autoclavable | | | | | | |
| 70 (2) | 2 3/4" | 3/8" | 2 11/16" | 12 | GY-06122-10 | |
| 125 (4) | 3 1/4" | 5/16" | 3 3/4" | 12 | GY-06122-20 | |
| 250 (8) | 4 1/4" | 7/16" | 4 11/16" | 12 | GY-06122-30 | |
| High-density polyethylene (HDPE) | | | | | | |
| 500 (16) | 5" | 1/2" | 5 1/8" | 12 | GY-06122-40 | |
| 1000 (32) | 6 1/2" | 15/16" | 7 1/8" | 12 | GY-06122-50 | |
| 2000 (64) | 8 1/8" | 3/4" | 8" | 12 | GY-06122-60 | |

Cole-Parmer® Funnels, 304 SS (medium gauge)

No wimpy construction—stands up to your needs

- Tab for easy lifting



| Capacity mL (oz) | Top OD | Total height | Catalog number | Qty/pk | Price/pk |
|------------------|--------|--------------|-----------------------------|--------|----------|
| 70 (2 3/8) | 3" | 3 1/4" | GY-07263-00 | 1 | |
| 190 (6 3/8) | 4 1/8" | 4 3/8" | GY-07264-00 | | |
| 400 (13) | 5" | 5 3/8" | GY-07265-00 | | |
| 650 (22 1/4) | 6" | 6 3/8" | GY-07266-00 | 1 | |
| 1000 (32) | 6 3/4" | 7 1/4" | GY-07267-00 | | |
| 2000 (64) | 8 3/8" | 8 3/8" | GY-07268-00 | | |

Brand 6220 Short-Stemmed Funnels

Easy sample transfer

- Clean up easily
- Ideal for powder transfer or bottle filling



Filtering Funnels, PP

No air lock
with outer ribs

- Get faster filtration with 58° angled internal ribs
- Accept standard filter papers
- Autoclavable



| Capacity mL (oz) | Approximate dimensions (mm) | | Catalog number | Qty | Price/ qty |
|---------------------|-----------------------------|----------------|-----------------------------|-------|---------------|
| | Top inside dia | Overall length | | | |
| 75 (2½) | 65 | 83 | GY-34564-01 | 12/pk | |
| 113 (3¾) | 75 | 94 | GY-34564-02 | 12/pk | |
| 261 (8¾) | 100 | 120 | GY-34564-03 | 6/pk | |
| 437 (14¾) | 125 | 132 | GY-34564-04 | 6/pk | |
| 748 (25¼) | 150 | 157 | GY-34564-05 | 6/pk | |
| 75 (2½) | 65 | 83 | GY-34564-21 | 24/cs | |
| 113 (3¾) | 75 | 94 | GY-34564-22 | 24/cs | |
| 261 (8¾) | 100 | 120 | GY-34564-23 | 24/cs | |
| 437 (14¾) | 125 | 132 | GY-34564-24 | 12/cs | |
| 748 (25¼) | 150 | 157 | GY-34564-25 | 12/cs | |

| Capacity mL (oz) | Top OD | Stem OD | Paper dia | Catalog number | Qty/ pk | Price/ pk |
|---------------------|-----------|------------|--------------|-----------------------------|------------|--------------|
| 10 (⅓) | 41 mm | 6 mm | 55 mm | GY-06120-10 | 12 | |
| 23 (¾) | 53 mm | 7 mm | 70 mm | GY-06120-15 | 12 | |
| 41 (1½) | 61 mm | 7 mm | 90 mm | GY-06120-20 | 12 | |
| 65 (2) | 72 mm | 7 mm | 110 mm | GY-06120-30 | 12 | |
| 114 (3¾) | 84 mm | 7 mm | 125 mm | GY-06120-40 | 6 | |
| 173 (6) | 98 mm | 8 mm | 150 mm | GY-06120-50 | 6 | |
| 254 (8) | 108 mm | 9 mm | 185 mm | GY-06120-60 | 4 | |
| 976 (32) | 166 mm | 14 mm | 240 mm | GY-06120-70 | 2 | |

Brand 6140 Filtering Funnels

Better accuracy with
filter paper

- 60° plain angle, long stem
- Beaded edge for additional strength
- Stems are manufactured with heavy-wall tubing



Brand 6180 Filtering Funnels

Filter faster
with depressed
internal flutings

- Feature fluted, 60° angle



| Capacity mL (oz) | Approximate dimensions (mm) | | Catalog number | Qty | Price/ qty |
|---------------------|-----------------------------|-------------------|-----------------------------|-------|---------------|
| | Top inside dia | Overall length | | | |
| 28 (1) | 50 mm | 193 mm | GY-34560-50 | 6/pk | |
| 62 (2) | 65 mm | 206 mm | GY-34560-51 | 6/pk | |
| 94 (3⅓) | 75 mm | 214 mm | GY-34560-52 | 6/pk | |
| 225 (7¾) | 100 mm | 236 mm | GY-34560-53 | 3/pk | |
| 28 (1) | 50 mm | 193 mm | GY-34560-60 | 12/cs | |
| 62 (2) | 65 mm | 206 mm | GY-34560-61 | 12/cs | |
| 94 (3⅓) | 75 mm | 214 mm | GY-34560-62 | 12/cs | |
| 225 (7¾) | 100 mm | 236 mm | GY-34560-63 | 6/cs | |

| Capacity mL (oz) | Approximate dimensions (mm) | | Catalog number | Qty | Price/ qty |
|---------------------|-----------------------------|----------------|-----------------------------|-------|---------------|
| | Top inside dia | Overall height | | | |
| 28 (1) | 50 | 108 | GY-34563-01 | 12/pk | |
| 62 (2) | 65 | 121 | GY-34563-02 | 12/pk | |
| 94 (3⅓) | 75 | 139 | GY-34563-03 | 12/pk | |
| 235 (8) | 100 | 190 | GY-34563-04 | 6/pk | |
| 28 (1) | 50 | 108 | GY-34563-21 | 48/cs | |
| 62 (2) | 65 | 121 | GY-34563-22 | 48/cs | |
| 94 (3⅓) | 75 | 139 | GY-34563-23 | 48/cs | |
| 235 (8) | 100 | 190 | GY-34563-24 | 24/cs | |

Cole-Parmer® Analytical Funnels, PTFE

Purity assured—chemically inert and nonstick

- Withstand most chemical compounds and corrosive environments
- Temperature usage is 32 to 536°F (0 to 280°C)
- Autoclavable



| Capacity mL (oz) | Top OD | Stem OD | Total height | Catalog number | Qty/ pk | Price/ pk |
|---------------------|-----------|------------|-----------------|-----------------------------|------------|--------------|
| 5 (⅓) | 30 mm | 8 mm | 50 mm | GY-06300-55 | 1 | |
| 23 (¾) | 50 mm | 10 mm | 85 mm | GY-06300-56 | 1 | |
| 70 (2) | 70 mm | 12 mm | 110 mm | GY-06300-57 | 1 | |
| 195 (7) | 100 mm | 14 mm | 170 mm | GY-06300-58 | 1 | |
| 700 (22) | 150 mm | 16 mm | 220 mm | GY-06300-59 | 1 | |



Funnels

Cole-Parmer® Plastic Büchner Funnels, HDPE and PP

Get quick, economical vacuum filtration

- Easy sterilization: HDPE by gas, PP by autoclave or gas



| Filter size (mm) | Catalog number | Qty/pk | Price/pk |
|---|-----------------------------|--------|----------|
| High-density polyethylene (HDPE) | | | |
| 42.5 | GY-06128-02 | 1 | |
| 55 | GY-06128-03 | 1 | |
| Polypropylene (PP), autoclavable | | | |
| 90 | GY-06128-04 | 1 | |
| 110 | GY-06128-05 | 1 | |
| 150 | GY-06128-10 | 1 | |
| 240 | GY-06128-20 | 1 | |

Cole-Parmer® Büchner Funnels, Porcelain

Make filtering easy

- Glazed with exception of rim
- Acid and alkali resistant
- Autoclavable




| Capacity mL (oz) | Filter dia (mm) | Height (mm) | Catalog number | Qty/pk | Price/pk |
|------------------|-----------------|-------------|-----------------------------|--------|----------|
| 50 (1¾) | 50 | 65 | GY-17800-00 | 1 | |
| 120 (4) | 70 | 76 | GY-17800-02 | 1 | |
| 200 (7) | 90 | 102 | GY-17800-04 | 1 | |
| 400 (13) | 110 | 120 | GY-63100-06 | 1 | |

Cole-Parmer® Separatory Funnels, PP

Won't leak when inverted

- Save money versus glass
- Leakproof threaded screw cap keeps a tight fit
- Self-lubricating PTFE stopcock plug



| Capacity mL (oz) | Closure size | Stem length | Stem OD | Catalog number | Qty/pk | Price/pk |
|------------------|--------------|-------------|---------|-----------------------------|--------|----------|
| 100 (3) | 28 mm | 66 mm | 8.5 mm | GY-06124-02 | 2 | |
| 250 (8) | | | | GY-06124-04 | | |
| 500 (16) | | | | GY-06124-06 | | |

Büchner Funnels, Ceramic

Stand up to harsh chemicals

- Glazed inside and out for easy cleaning
- Autoclavable



| Capacity (mL) | Filter dia (mm) | Catalog number | Qty/pk | Price/pk | Catalog number | Qty/cs | Price/cs |
|---------------|-----------------|-----------------------------|--------|----------|-----------------------------|--------|----------|
| 2 | 13 | GY-17815-02 | 1 | | GY-17815-22 | 12 | |
| 30 | 35–40 | GY-17815-04 | 1 | | GY-17815-24 | 12 | |
| 87 | 50–55 | GY-17815-06 | 1 | | GY-17815-26 | 8 | |
| 186 | 70 | GY-17815-08 | 1 | | GY-17815-28 | 6 | |
| 320 | 90 | GY-17815-10 | 1 | | GY-17815-30 | 4 | |
| 550 | 110 | GY-17815-12 | 1 | | GY-17815-32 | 4 | |
| 700 | 110–124 | GY-17815-14 | 1 | | GY-17815-34 | 2 | |
| 1060 | 142–150 | GY-17815-16 | 1 | | GY-17815-36 | 2 | |
| 1860 | 185 | GY-17815-18 | 1 | | GY-17815-38 | 2 | |
| 4500 | 240 | — | — | — | GY-17815-40 | 1 | |

Brand 6404 Separatory Funnels

Easy-to-grasp HDPE stopper

- Eliminate contamination—PTFE stopcock requires no grease
- Long 65-mm stem for easier sampling in deep containers
- Case quantities also available



| Capacity mL (oz) | § stopper no. | § stopper plug size | Approx OD (mm) | Approx height (mm) | Catalog number | Qty/pk | Price/pk |
|------------------|---------------|---------------------|----------------|--------------------|-----------------------------|--------|----------|
| 60 (2) | 16 | 2 | 10 | 244 | GY-34507-30 | 1 | |
| 125 (4) | 22 | 2 | 10 | 289 | GY-34507-32 | | |
| 250 (8) | 22 | 4 | 10 | 311 | GY-34507-34 | | |
| 500 (16) | 27 | 4 | 10 | 354 | GY-34507-36 | 1 | |
| 1000 (32) | 27 | 4 | 10 | 392 | GY-34507-58 | | |
| 2000 (64) | 38 | 6 | 12 | 486 | GY-34507-60 | | |

§ Symbol used to designate interchangeable joints, stoppers and stopcocks that comply with the requirements of Commercial Standard CCU-21 published by NIST.

§ Symbol designates stopcock plugs made of PTFE that meet requirements of NIST Voluntary Product Standard PS 28-70.

Wait!

There's More at ColeParmer.com

Find more than **100,000** products to choose from!



Plastic Separatory Funnels, PP and FEP

No leaks or lubrication

- Screw cap keeps contents inside funnel

| Capacity mL (oz) | Closure size | Stem length | Stem OD | Catalog number | Qty/pk | Price/pk |
|---|--------------|-------------|---------|-----------------------------|--------|----------|
| Polypropylene (PP) body with polypropylene screw closure | | | | | | |
| 125 (4) | 28 mm | 65 mm | 8.5 mm | GY-06126-10 | 1 | |
| 250 (8) | 33 mm | | | GY-06126-20 | | |
| 500 (16) | 43 mm | | | GY-06126-50 | | |
| 1000 (32) | 53 mm | | | GY-06126-60 | | |
| FEP body with ETFE screw closure | | | | | | |
| 125 (4) | 28 mm | 65 mm | 8.5 mm | GY-06125-10 | 1 | |
| 250 (8) | 33 mm | | | GY-06125-20 | | |
| 500 (16) | 43 mm | | | GY-06125-50 | | |
| 1000 (32) | 53 mm | | | GY-06125-60 | | |
| 2000 (64) | 53 mm | | | GY-06125-70 | | |



Separatory Funnels

Won't leak when shaken or transported

- Externally threaded outer ground joint secures stopper with a cap and O-ring
- Include PTFE plug, stopper, cap, and O-ring

| Capacity, mL (oz) | Stopper size | Stopcock no. | Catalog number | Qty/cs | Price/cs |
|-------------------|--------------|--------------|-----------------------------|--------|----------|
| 60 (2) | 24/40 | 2 | GY-34731-00 | 1 | |
| 125 (4) | 24/40 | | GY-34731-02 | | |
| 250 (8) | 24/40 | 4 | GY-34731-04 | 1 | |
| 500 (16) | 24/40 | | GY-34731-06 | | |
| 1000 (32) | 24/40 | 4 | GY-34731-08 | 1 | |
| 2000 (64) | 45/50 | | GY-34731-10 | | |



Brand 6402 Separatory Funnels

Excellent chemical resistance

- Pyrex stopper, PTFE stopcock
- No freezing or grease contamination from stopcock
- Case quantities also available



Antistatic Plastic Weighing Funnels, PP

No more messy transfers to volumetric flasks

- No-roll, flat side for placement on balance
- Stem fits easily in volumetric flasks for easier transfer



Flat-bottom funnels are stable and easy to weigh.

| Capacity mL (oz) | Stopper no. | Stopper plug size | Approx height (mm) | Catalog number | Qty/pk | Price/pk |
|------------------|-------------|-------------------|--------------------|-----------------------------|--------|----------|
| 60 (2) | 16 | 2 | 225 | GY-34506-01 | 1 | |
| 125 (4) | 22 | 2 | 280 | GY-34506-02 | | |
| 250 (8) | 22 | 4 | 311 | GY-34506-03 | | |
| 500 (16) | 27 | 4 | 355 | GY-34506-04 | | |
| 1000 (32) | 27 | 4 | 392 | GY-34506-25 | | |
| 2000 (64) | 38 | 6 | 496 | GY-34506-26 | 1 | |
| 4000 (128) | 38 | 8 | 575 | GY-34506-27 | | |
| 6000 (200) | 38 | 10 | 600 | GY-34506-28 | | |

| Size | Fits flask size (mL) | Neck OD | Cat. no. | Qty/pk | Price/pk |
|-------------|----------------------|---------|-----------------------------|--------|----------|
| Micro | 0.1 to 1† | 6 mm | GY-06127-00 | 96 | |
| Small | 5, 10, 25, and 50 | 7 mm | GY-06127-02 | 50 | |
| Medium | 100, 200, and 250 | 11 mm | GY-06127-04 | 50 | |
| Large | 500 and 1000 | 18 mm | GY-06127-06 | 20 | |
| Extra large | 1000 and 2000 | 20 mm | GY-06127-08 | 10 | |

†Micro funnel fits 96-well mini-tubes and deep well plates.

Separatory Funnels

Get a more complete separation

- Larger stem ID allows separation even when PTFE stopcock is closed
- Meets ASTM specification E1096, Type 4 requirements



| Capacity mL (oz) | Stopper no. | Stopper plug size | Stem length | Catalog number | Qty/cs | Price/cs |
|------------------|-------------|-------------------|-------------|-----------------------------|--------|----------|
| 30 (1) | 13 | 2 | 60 mm | GY-34505-40 | 4 | |
| 60 (2) | 16 | 2 | 60 mm | GY-34505-42 | 4 | |
| 125 (4) | 22 | 2 | 60 mm | GY-34505-44 | 4 | |
| 250 (8) | 22 | 4 | 60 mm | GY-34505-46 | 4 | |
| 500 (16) | 27 | 4 | 60 mm | GY-34505-48 | 4 | |
| 1000 (32) | 27 | 4 | 60 mm | GY-34505-50 | 2 | |
| 2000 (64) | 38 | 6 | 60 mm | GY-34505-52 | 2 | |

Glass Weighing Funnels

Save money—reusable

- Ensure complete sample transfer to volumetric flasks
- Flat on one side for accurate weighing on a balance
- Micro funnel fits microwell plates



Flat-bottom design will not roll.

| Size | Fits flask size (mL) | Neck OD | Catalog number | Qty/pk | Price/pk |
|-------------|----------------------|---------|-----------------------------|--------|----------|
| Micro | 0.1 to 1† | 7 mm | GY-34700-00 | 6 | |
| Small | 5, 10, 25, 50 | 8 mm | GY-34700-02 | 6 | |
| Medium | 100, 200, 250 | 11 mm | GY-34700-04 | 6 | |
| Large | 500, 1000 | 17 mm | GY-34700-06 | 6 | |
| Extra large | 1000, 2000 | 20 mm | GY-34700-08 | 6 | |

†Micro funnel fits 96-well mini-tubes and deep well plates.