


Racks

Economical Storage Boxes, PP

Won't deteriorate in moisture

- Use in temperatures down to -80°C (-112°F)
- Easily identify samples through translucent lids
- Autoclavable



Tube size (mL)	Tubes held	Maximum tube length (mm)	Qty/pk	Catalog number	Price/pk
Boxes with unattached lid					
0.5	50	31	6	GY-44350-00	
1.5	50	43	6	GY-44350-02	
Box with hinged lid					
1.5	100	49	6	GY-44350-10	

Switch-Grid™ Tube Racks

Get 25 racks in one!

- One rack provides up to 25 possible tube size combinations with quick changes
- Accommodates up to two different-sized tubes in the same rack
- Less expensive than purchasing multiple racks
- Chemical-resistant design provides easy cleaning
- Autoclavable



Tube size (mm)	No. of tubes held (array)	Catalog number	Price
10 to 13	84 (14 x 6)	GY-06736-12	
13 to 16	60 (12 x 5)	GY-06736-13	
16 to 20	40 (10 x 4)	GY-06736-14	
20 to 25	24 (8 x 3)	GY-06736-16	
25 to 30	12 (6 x 2)	GY-06736-17	

Accessories

- [GY-06736-18](#) Grid inserts for 10 to 13 mm tubes; orange. Set of 4
- [GY-06736-19](#) Grid inserts for 13 to 16 mm tubes; blue. Set of 4
- [GY-06736-22](#) Grid inserts for 16 to 20 mm tubes; green. Set of 4
- [GY-06736-23](#) Grid inserts for 20 to 25 mm tubes; fuchsia. Set of 4
- [GY-06736-24](#) Grid inserts for 25 to 30 mm tubes; yellow. Set of 4

Test Tube Support System, PP

One rack, virtually limitless versatility

- Rotate cubes to accommodate different size tubes
- Consists of four cubes and tray
- Autoclavable



Compatible tube sizes	Catalog number	Price
15- and 50-mL conical; 12 x 75 mm, 13 x 100 mm, 16 x 125 mm, and 17 x 100 mm culture; and 1.5- and 2.0-mL microcentrifuge	GY-06733-00	

Test Tube Racks, Epoxy-coated steel

Won't rust or corrode

- Epoxy-coated steel wire racks
- Unaffected by organic solvents, salts, weak acids and alkalies
- Handle temperatures from -369 to 250°F (-223 to 121°C)
- Autoclavable



Tube size (mm)	Tubes held (array)	Qty/pk	Catalog number	Price/pk
10 to 13	72 (6 x 12)	1	GY-06739-02	
15 to 16	48 (4 x 12)	1	GY-06739-04	
18 to 20	40 (4 x 10)	1	GY-06739-06	
22 to 25	40 (4 x 10)	1	GY-06739-08	

Lab-Aire® II Drying Racks

Configure to your labware

- Reposition pegs to accommodate labware
- Scalloped well drip tray allows for drying of burettes and pipettes
- Won't corrode

Description	No. of pegs	Dimensions (W x H x D)		Catalog number	Price
Single-sided drying racks					
1 tier, benchtop/wall mount	19	16" x 9¼" x 7½"	40.6 x 23.5 x 19.1 cm	GY-06733-91	
2 tier, benchtop	38	16" x 16¼" x 7½"	40.6 x 41.3 x 19.1 cm	GY-06733-93	
2 tier, wall mount				GY-06733-98	
3 tier, benchtop	57	16" x 23¾" x 7½"	40.6 x 60.3 x 19.1 cm	GY-06733-94	
3 tier, wall mount				GY-06733-99	
Double-sided drying racks					
1 tier, benchtop	38	16" x 9¼" x 10½"	40.6 x 23.5 x 26.7 cm	GY-06733-95	
2 tier, benchtop	76	16" x 16¼" x 10½"	40.6 x 41.3 x 26.7 cm	GY-06733-96	
3 tier, benchtop	114	16" x 23¾" x 10½"	40.6 x 60.3 x 26.7 cm	GY-06733-97	




Wall-mount racks are sized to fit standard pegboards

Cole-Parmer® Drying/Draining Rack, PP

Don't let water collect

- Draining port allows water to drain back to sink
- Pegs screw-in for removal to accommodate different size labware
- Features 5 tiers; 20 pegs



Dimensions (W x H x D)		Catalog number	Price
13¾" x 19¼" x 2¼"	34.9 x 48.9 x 5.7 cm	GY-06045-91	

Cole-Parmer® Drying Rack, PS

Dry more glassware

- Includes 72 v-shaped pegs to accommodate more labware
- Easy snap in peg design



Rack 06045-95 shown with optional acrylic stand 06045-98



Dimensions (W x H x D)		Catalog number	Price
17¾" x 24" x 3"	45.1 x 61 x 7.6 cm	GY-06045-95	

[GY-06045-97](#) Optional test tube pegs, rounded. Bag of 11
[GY-06045-98](#) Optional acrylic stand with mounting hardware

Test Tube Racks, PP

Hold samples in baths effectively

- Open rack drains fluids quickly
- Models 06737-10, -15, and -20 stack with ease
- Autoclavable



Tube size (mm)	Tubes held (array)	Qty/pk	Catalog number	Price/pk
10 to 13	72	1	GY-06737-10	
13 to 16	40	1	GY-06737-15	
16 to 20	40	1	GY-06737-20	
21 to 25	40	1	GY-06737-30	
26 to 30	24	1	GY-06737-40	

Labware Drain Stand, Epoxy-Coated Steel/PS

Epoxy-Coated Steel/PS

Dry various size labware

- Accommodates up to 26 pieces of labware with two different size pins
- No more mess with convenient drain tray



Dimensions (W x H x D)		Catalog number	Price
16½" x 18½" x 6½"	43.0 x 48.0 x 17.0 cm	GY-06739-99	

Modular Drying Racks, 304 SS

Don't settle—customize to your needs in seconds

- Easy-to-clean, corrosion-resistant stainless steel racks and polypropylene (PP) pegs
- No mess—catch drain directs drainage to drip trough
- Autoclavable
- Include 6" L pegs, SS drain tube, and wall-mounting bracket



Rack 67035-02 shown with drain basket and glove dispenser



Rack 67035-10 shown with pipette rack and screen insert




Color of pegs (included)	Catalog number	Price/ea
32-peg racks measure 24"W x 30"H x 4"D (61 x 76.2 x 10.2 cm)		
White	GY-67035-00	
Assorted†	GY-67035-02	
50-peg racks measure 30"W x 30"H x 4"D (76.2 x 76.2 x 10.2 cm)		
White	GY-67035-10	
Assorted†	GY-67035-12	
Black	GY-67035-14	
60-peg racks measure 36"W x 30"H x 4"D (91.4 x 76.2 x 10.2 cm)		
White	GY-67035-20	
Assorted†	GY-67035-22	
Black	GY-67035-24	

†Assorted peg colors include yellow, orange, pink, red, green, light blue, dark blue, lilac, and purple.

Add-On Accessories, 304 SS

These accessories include the necessary "snap-in" support pegs to quickly and easily attach to the rack. **Drain basket** holds small items such as stir bars, stoppers, and pipette tips. **Drain shelf** is ideal for drying and storing various glassware items. **Pipette rack** hangs at an angle to allow for complete drainage of pipettes. **Funnel rack** holds four funnels, volumetric flasks, or small cylinders. **Cylinder yoke** holds two cylinders. **Flask holder** accepts one large flask, cylinder, or separatory funnel. **Glove dispenser** provides a convenient spot for one box of gloves. **Screen inserts** fit into the drip trough and increase the drying capacity of the rack.



Drain basket 67035-70



Drain shelf 67035-71




Pipette rack 67035-72



Funnel rack 67035-73




Cylinder yoke 67035-74



Flask holder 67035-75



Glove dispenser 67035-76



Screen insert 67035-77

Description	Catalog number	Price/ea
Drain basket, 6"W x 2 1/2"H (15.2 x 6.4 cm)	GY-67035-70	
Drain shelf, 14"W (35.6 cm)	GY-67035-71	
Pipette rack, 12 1/2"L (31.8 cm)	GY-67035-72	
Funnel rack, 14"W (35.6 cm)	GY-67035-73	
Cylinder yoke, 14"W (35.6 cm)	GY-67035-74	
Flask holder, 8"W (20.3 cm)	GY-67035-75	
Glove dispenser, 10"W x 5 1/4"H (25.4 x 14.6 cm)	GY-67035-76	
Screen insert, 24"W (61 cm)	GY-67035-77	
Screen insert, 30"W (76.2 cm)	GY-67035-78	
Screen insert, 36"W (91.4 cm)	GY-67035-80	

Plugs and Pegs, Molded Polypropylene

Removable "snap-in" plugs and pegs come in packs of six.

Description	Catalog number	Price/pk
White hole plugs	GY-67035-79	
White tube pegs, 4"L x 5/16" dia (10.2 x 0.8 cm)	GY-67035-90	
White long pegs, 9"L x 5/8" dia (22.9 x 1.6 cm)	GY-67035-91	
White long pegs, 12"L x 5/8" dia (30.5 x 1.6 cm)	GY-67035-92	