

Tongs & Tweezers

Straight Blade Scissors, Stainless Steel and Ceramic

Cut through most materials

Choose the all-purpose stainless steel scissors for a variety of uses.


The ceramic straight blade is ideal for precision cuts and works with magnetic materials; features plastic and carbon fiber handles.

Description	Overall length	Catalog number	Qty/pk	Price/pk
Stainless steel	4½"	GY-07292-00	1	
	8"	GY-07292-10	1	
Ceramic	5"	GY-07292-20	1	



Tongs

Get a grip on your samples



Key	Description	Material	Length	Catalog number	Price/ea
A	Individually wrapped, sterile sampling tongs, serrated ends to securely hold wet and dry samples	Polystyrene	4¼" (10.8 cm)	GY-06443-02	
B	Tongs feature ¾" square tips for easy gripping; autoclavable	ETFE	10" (25.4 cm)	GY-06442-00	
C	Beaker tongs feature flat riveted joint, jaw opens to 7"; handles beakers from 50 to 1500 mL	Zinc-plated steel	9" (22.9 cm)	GY-06443-06	
D	Crucible tongs, 3" PTFE coating on each tip	Stainless steel	9" (22.9 cm)	GY-06443-00	
E	Tong is entirely covered in PTFE coating, no loss of sample through sticking; will not scratch or damage delicate surfaces	Stainless steel	9¾" (23.9 cm)	GY-06443-04	
F	Crucible tongs have positive-grip corrugated tips for safely handling dishes in muffle furnaces; holds crucibles or dishes 1½" or wider	Nickel-plated steel	20" (50.8 cm)	GY-06443-50	
G	Crucible tongs, riveted joint and serrated tips	Stainless steel	9" (22.9 cm)	GY-06443-07	
H	Beaker tongs with silicone tips, feature riveted joint; handles beakers up to 2 L	Stainless steel	12" (30.5 cm)	GY-06443-08	

Precision Tweezers, Titanium

Won't shed particles

Titanium tweezers are excellent for the silicon chip industry, for handling microscopic particles, and when working in a magnetic field as they are nonmagnetic, lightweight, and strong. They do not shed particles like stainless steel tweezers and are biocompatible. Maximum temperature resistance is 800°F (426°C).

Key	Description	Length	Catalog number	Price/ea
A	Titanium, style AA—Very strong, straight, precise tips	4¾" (12.1 cm)	GY-07289-20	
B	Titanium, style OO—Heavy duty strong blades that taper to regular strong points	4¾" (12.1 cm)	GY-07289-22	
C	Titanium, style 2A—Straight flat tips provide large flat gripping surface	4¾" (12.1 cm)	GY-07289-26	
D	Titanium, style 3C—Straight ultra fine tips	4¾" (12.1 cm)	GY-07289-30	
E	Titanium, style 7—Very sharp points	4¼" (10.8 cm)	GY-07289-36	



General-Purpose Tweezers and Forceps, Stainless Steel (SS) and Plastic

Get a more precise grip

SS Tweezers,
straight tips and finger serrations.



Length	Catalog number	Qty/ pk	Price/pk
7"	GY-07288-04	12	

SS Straight: Thick, Flat Strong Tips,
serrated tips and grips. Style OOD.



Length	Catalog number	Qty/ pk	Price/pk
4¾"	GY-07287-03	1	

SS Self-Locking Tweezers,
straight narrow tips.



Length	Catalog number	Qty/ pk	Price/pk
4¼"	GY-07288-01	12	

SS Tweezers,
smooth spade tips.



Length	Catalog number	Qty/ pk	Price/pk
4¼"	GY-07288-07	12	

SS Solder Utility Tweezers,
curved tips and finger serrations.



Length	Catalog number	Qty/ pk	Price/pk
7"	GY-07288-03	12	

SS Tweezers,
curved serrated tips.



Length	Catalog number	Qty/ pk	Price/pk
4¼"	GY-07288-05	12	

SS Forceps,
straight serrated tips.



Length	Catalog number	Qty/ pk	Price/pk
6"	GY-07288-12	1	
12"	GY-07288-14	1	

SS Forceps,
plastic-coated rounded tips.



Length	Catalog number	Qty/ pk	Price/pk
12"	GY-07288-10	1	

PTFE-Covered Nickel-Plated Forged Steel Forceps are heat and chemical resistant. **Autoclavable.**




Length	Catalog number	Qty/ pk	Price/pk
4¼"	GY-06440-00	2	

Reinforced Nylon Forceps have corrugated handles, serrated blades, and rounded tips. **Autoclavable.** Sterile—individually wrapped.




Sterile	Length	Catalog number	Qty/ pk	Price/pk
No	4⅞"	GY-06443-15	100	
Yes	4⅞"	GY-06443-26	100	

Nylon Forceps feature corrugated handles and serrated blades. **Autoclavable.** Sterile—individually wrapped.




Sterile	Length	Catalog number	Qty/ pk	Price/pk
No	4¼"	GY-06443-25	100	
Yes	4¼"	GY-06443-27	70	

PA (Polyamide) Tweezers,
rounded tips.



Length	Catalog number	Qty/ pk	Price/pk
7"	GY-06446-20	1	

Polystyrene Disposable Sterile Forceps feature corrugated handles and rounded serrated tips. Sterile—come individually wrapped.



Length	Catalog number	Qty/ pk	Price/pk
4⅞"	GY-06443-20	100	

Plastic Curved: Smooth Rounded Tips 2AB



Length	Catalog number	Qty/ pk	Price/pk
4½"	GY-07290-06	12	

Plastic Curved: Very Fine and Smooth Tips 7-STD



Length	Catalog number	Qty/ pk	Price/pk
4½"	GY-07290-08	12	